

References

1. Bach G, Neckel C, Mall C, Krekeler G. Conventional Versus Laser-Assisted Therapy of Periimplantitis: A Five-Year Comparative Study. *Implant Dent* 2000;9:247-251 **Considerably effective in vivo**
2. Romanos G, Javed F, Delgado-Ruiz R, Calvo-Guirado J. Peri-implant Diseases A Review of Treatment Interventions. *Dental Clinics of North America*:2015, 59 (1), 157-178 **Effective in vivo, temperature increase a major risk with CO2 and Er:YAG after 10 seconds. Nd:YAG laser is contraindicated**
3. Bach G, Diode laser surface decontamination in periodontitis therapy 15 years of incorporating. *Int M of Oral Implantology* 2011; 3 (1) 7-12 **Considerably effective, no thermal damage in vivo**
4. Bach G, Neckel C, Frejbu G. A 5-year comparative study on conventional and laser assisted therapy of periimplantitis and periodontitis. *Implant Dent*. 2000 9, 247-251 **Significantly effective**
5. Romanos G, Nentwig G, Regenerative Therapy of Deep Peri-implant Infrabony Defects After CO2 Laser Implant Surface Decontamination. *Int J Periodontics Restorative Dent* 2008;28:245-255. **Considerably effective, diode only alternative to CO2**
6. Pai J, Shridhar A, Kamath V, Jaiswal N, Malagi S, Radhika B. Adjunctive use of diode lasers in the treatment of peri-implantitis: A case series. *Journal of Dental Implants*. 2014;4(1):91-97 **Safe effective positive biostimulation effects in vivo in contact and non contact**
7. Mettraux G, A concept of laser assisted treatment of periimplantitis *Int M of Oral Implantology* 2011; 3 (1) 30-35 **Obvious bone regeneration after 2 years, biofilm management**
8. Roncati M, Lucchese A, Carinci F. Non-surgical treatment of peri-implantitis with the adjunctive use of an 810-nm diode laser. *JIS of Periodontology*.2013;17(6):812. **Effective non surgical**
9. Cobb CM. Lasers in periodontics: A review of the literature. *J Periodontal*. 2006;77:545-64. **Diode lasers have a bactericidal effect, 278 possible peer reviewed articles.**
10. Moritz A, Schoop U, Goharkay K, Shauer P, Doertbuduk O, Wernisch J, et al. Treatment of periodontal pockets with a diode laser. *Lasers Surg Med* 1998;22:302-11. **Bacterial reduction with diode laser therapy was significantly better**
12. Romanos GE, Treatment of periimplant lesions using different laser systems. *Oral Laser Appl*. 2002; 2; 75-80 **Superficial absorption with no damage to underlying tissue with diode**
13. Romanos G, Nentwig GH, Diode Laser (980 nm) in Oral and Maxillofacial Surgical Procedures: Clinical Observations Based on Clinical Applications. *Journal of Clinical Laser Medicine and Surgery*. 1999; 17 (5) 193-197 **Precise incision margin**
14. Gonçalves F¹, Zanetti AL, Zanetti RV, Martelli FS, Avila-Campos MJ, Tomazinho LF, Granjeiro JM. Effectiveness of 980-nm diode and 1064-nm extra-long-pulse neodymium-doped yttrium aluminum garnet lasers in implant disinfection. *Photomed Laser Surg*. 2010; 28(2):273-80. **100% bacteria reduction without damaging implant**
15. Tosun E, Tasar F, Strauss R, Kivanc G, Ungor C. Comparative Evaluation of Antimicrobial Effects of Er:YAG, Diode, and CO2 Lasers on Titanium Discs: An Experimental Study *J Oral Maxillofac Surg* 2012; 70:1064-1069 **100% bacteria elimination with diode**
16. Romanos G, Hayo G, Nentwig G. Effects of Diode and Nd:YAG Laser Irradiation on Titanium Discs: A Scanning Electron Microscope Examination. *J Periodontol* 2000;71:810-815 **980 does not damage implant regardless of power setting**
17. Garg H, Garg V, Jagadeesh H G, Bedi G. Peri-Implantitis and Lasers-A Review. *Int J Biol Med Res*. 2012; 3(3):2302-2305 **Significant reduction bacteria in pockets, no damage with 980nm diode. Does not recommend the use of CO2 or Nd:YAG lasers in peri implantitis.**
18. Chris Leja C, Geminiani A, Caton J, Romanos G. Thermodynamic effects of laser irradiation of implants placed in bone: an in vitro study. *Lasers Med Sci*. 2013 Nov;28(6):1435-40. **980 pulsed was only laser to not reach critical threshold**
19. Khandge N, Pradhan S, Doshi Y, Kulkarni A. Comparison of the Effects of Different Laser Wavelengths on Implants Surfaces. *International Journal of Laser Dentistry*, 2013;3(1):14-1 **980 pulsed most effective and only laser without negative surface defects, better than 810 performance, coagulation and deeper incisions than CO2 or Nd:YAG laser at the same power setting**
20. Hauser-Gerspach I¹, Stübinger S, Meyer J. Bactericidal effects of different laser systems on bacteria adhered to dental implant surfaces: an in vitro study comparing zirconia with titanium. *Clin Oral Implants Res*. 2010 Mar;21(3):277-83. **Effectively reduced bacteria**
21. Goharkhay K, Moritz A, Wilder-Smith P, Schoop U, Kluger W, Jakolitsch S, Sperr W. Effects on oral soft tissue produced by a diode laser in vitro. *Lasers Surg Med* 1999; 25:401-406 **Pulsed eliminates risk of thermal damage, safer than CW.**
23. Kreisler M¹, Götz H, Duschner H. Effect of Nd:YAG, Ho:YAG, Er:YAG, CO2, and GaAlAs (Diode) laser irradiation on surface properties of endosseous dental implants. *Int J Oral Maxillofac Implants*. 2002;17(2):202-11. **Nd:YAG and Ho:YAG lasers are not suitable. Er:YAG and CO2 limit power to avoid damage. Diode safe.**
25. Romanos GE, Crespi R, Barone A, Covani U. Osteoblast attachment on titanium disks after laser irradiation. *Int J Oral Maxillofac Implants* 2006;21:232-236. **Osteoblasts can be grown, may promote osteoblast attachment**
26. Walsh LJ. The use of lasers in implantology: An overview. *J Oral Implantol* 1992;18:335-40. **Diode effective, Nd:YAG not recommended**
30. Romanos GE, Treatment of periimplant lesions using different laser systems. *Oral Laser Appl*. 2002; 2; 75-80 **No damage to underlying tissue**
31. Valderrama P, Blansett J, Gonzalez M, Cantu M, Wilson T. Detoxification of Implant Surfaces Affected by Peri-Implant Disease: An Overview of Non-surgical Methods. *The Open Dentistry Journal*, 2014,8, (Suppl 1-M5) 77-84 **Co2 induced surface changes in the hydroxyapatite. Er:YAG laser caused surface alterations. Diode sterilized implant surfaces. Nd:YAG and Ho:YAG lasers not recommended**
32. Oppenheimer A. Temperature Changes Resulting from a GaAlAs (Diode) Laser in the Decontamination of a Failing Dental Implant. *Jour of S Car Acad of Sci*. 2009; (7) 2 28 **CW exceeds temperature threshold. Pulse mode allows for higher power setting without reaching threshold for optimal bacteria kill**
34. Bains K, Gupka S, Bains R, Lasers in Periodontics: An Overview. *JOHCD*.2010; 4 29 **Nd:YAG laser is contraindicated, CO2 risk with high temperature**
35. Park CY¹, Kim SG, Kim MD, Eom TG, Yoon JH, Ahn SG. Surface properties of endosseous dental implants after NdYAG and CO2 laser treatment at various energies. *J Oral Maxillofac Surg*. 2005; 63(10): 1522-7. **Nd:YAG not recommended**
36. Yousif A, Zwinger S, Beer F, Verhagen L, Iwitschier M, Strassle M, Wintner E. Investigation on Laser Dental Implant Decontamination. *JLMN-Journal of Laser Micro/Nanoengineering*. 2008; Vol. 3, No. 2, **Diode must be calibrated. Diode lasers real output parameters can differ significantly.**
37. Romanos G, Nentwig GH, Diode Laser (980 nm) in Oral and Maxillofacial Surgical Procedures: Clinical Observations Based on Clinical Applications. *Journal of Clinical Laser Medicine and Surgery*. 1999; 17 (5) 193-197 **Precise incision margin, coagulation of vascular lesions, lack of swelling, bleeding, pain or, scar tissue formation, and good wound healing were observed in all of the clinical applications.**

18+ Years, Over 5000 Owners, 22 Residency Programs